



## LiteScope™ laryngoscope

One adult size ▪ Fully disposable ▪ Ready for use



**Airway Management** ▪ Airway Accessories

## Introducing the LiteScope™ laryngoscope

### One adult size • Fully disposable • Ready for use

LiteScope is the new, one piece, single-use, fully disposable laryngoscope from Intersurgical, offering a quick and easy-to-use solution for intubation whilst maximising patient safety.

Unlike a standard laryngoscope, the innovative one-piece design of LiteScope eliminates the need to connect a separate blade and handle, thus saving valuable time in an emergency. Modelled on a Macintosh size 4 blade, the insertion technique is the same as for any standard laryngoscope with a Macintosh blade.

LiteScope is ideal for use in:

- Pre-hospital
- Emergency medicine
- Resuscitation
- Anaesthesia
- Difficult airway
- Intensive care

### Infection Control

Recently published infection control guidelines have confirmed that in relation to laryngoscopes, 'current practices for decontamination and disinfection between patients are frequently ineffective, leaving residual contamination that has been implicated as a source of cross infection'. Blades are also regularly contaminated with blood, indicating penetration of mucous membranes, and puts these items into a high-risk category. They also confirm that, 'Ideally, a fully disposable laryngoscope, whereby the handle is also discarded with the blade, is recommended.'


As a single use, fully disposable laryngoscope, LiteScope eliminates the risks associated with the reprocessing of reusable devices, and therefore provides an ideal solution where there is a concern regarding infection control. It also avoids any reprocessing costs.



Anaesthesia



Intensive care



Emergency medicine

## Features & benefits

### One adult size

no need to stock multiple sizes

### Fully disposable

no need to remove batteries prior to disposal

### Ready for use

No need to connect separate blades

### Corrugated handle

For superior grip



### On/Off switch

LiteScope is ready to use seconds after removing from the packaging

### TPE guard

May help reduce the risk of dental trauma during intubation

### Macintosh blade

Familiar insertion technique

### Ergonomic design


Ensures LiteScope™ is comfortable and easy to use

### LED light

For optimal visualisation



[Watch the video](#)



LiteScope™ is supplied in boxes of 10 units. Each device is individually packed to ensure protection of the device, both in transit and in storage.



[www.intersurgical.com/info/litescope](http://www.intersurgical.com/info/litescope)

Code	Description	Box Qty.
7124000	LiteScope laryngoscope	10


[Make an enquiry](#)

**Intersurgical range of airway devices and accessories**

The LiteScope is the most recent addition to our range of airway devices and accessories, which includes the innovative i-gel® 2nd generation supraglottic airway, the i-view™ video laryngoscope and a whole range of laryngeal masks, laryngoscope blades, tracheal tubes, oropharyngeal airways and intubating bougies and stylets. It also includes an extensive range of breathing filters, heat and moisture exchangers (HMEs) and heat and moisture exchanging filters (HMEFs). For further details regarding any of our airway products, please do not hesitate to contact us.




i-gel® supraglottic airway information sheet




i-view™ video laryngoscope information sheet



InterGuide tracheal tube introducers and InterForm intubating stylets information sheet



Safe and effective breathing filters information sheet



One-piece Guedel airways information sheet

IS6.17 BENL • Issue 1 03.21



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland  
 T: +31 (0)413 243 860 F: +31 (0)413 243 869 [info@intersurgical.be](mailto:info@intersurgical.be) [www.intersurgical.be](http://www.intersurgical.be)



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016, ISO 14001:2015 and MDSAP

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.