



# nFlow™ CosyFit

Infant nasal CPAP system



**Critical Care** ▪ Breathing Systems



# nFlow™ CosyFit

## Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

### A choice for infant CPAP

The Intersurgical nFlow™ system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.

### The CosyFit bonnet

The CosyFit bonnet makes attaching the generator simple and precise and helps create a secure yet gentle fit for even the most sensitive patients.

### Tube support



Lightweight and secure design for easy positioning and adjustment

### Bonnet tie

Soft and strong to maintain a comfortable and reliable fit throughout treatment

### Bonnet

Available in 12 sizes, the breathable cotton nFlow CosyFit bonnet includes 3 large attachment pads to help ensure a gentle and cosy fit for the smallest to largest patient.



### The nFlow system


The nFlow system is supplied assembled with generator, heated wire limb and 10mm adjustable expandable exhaust tube for efficient gas removal. The tube features pressure safety slits to prevent the build up of patient pressure in the event of occlusion. The complete nFlow breathing system is phthalate free\* including the small bore tube used for monitoring and gas delivery.

\* excluding product code 4861000



## A choice of patient interfaces and bonnets

Patient interfaces and CosyFit bonnets are supplied separately to avoid waste. Use the size guide provided to choose the correct option for your patient.




### 12 bonnets

Available separately in twelve colour-coded sizes. Each CosyFit bonnet includes 3 fixation pads to secure the tube support and bonnet ties.



### Patient interface

Masks and prongs are available separately in four size options. Made from soft, comfortable silicone to ensure optimal patient fit with minimal wastage.



## B CosyFit bonnets



## C Nasal masks



## D Nasal prongs




Make an enquiry

Find out more

▶ Watch the video



We've made ordering easy. To make a complete system you can choose one product from group A and B, and choose a patient interface from C and/or D. At any time you can order separate items according to your needs.



Code	Description	Box Qty.	
<b>A nFlow CosyFit breathing systems</b>			
4800000	nFlow CosyFit infant CPAP heated wire breathing system. ≥ 1.6m	12	
4800001	nFlow CosyFit infant CPAP heated wire breathing system. ≥ 1.8m	10	
4800016	nFlow CosyFit infant CPAP heated wire breathing system with detachable monitoring line. ≥ 1.8m	10	
4861000	nFlow CosyFit infant CPAP heated wire breathing system for Fabian™ ventilators with detachable monitoring line. ≥ 1.8m	10	
<b>B nFlow CosyFit bonnets</b>			
		<b>Size</b>	
4907009	nFlow CosyFit, infant bonnet, navy	Size 9 (40 - 42 cm)	10
4907008	nFlow CosyFit, infant bonnet, turquoise	Size 8 (38 - 40 cm)	10
4907007	nFlow CosyFit, infant bonnet, orange	Size 7 (36 - 38 cm)	10
4907006	nFlow CosyFit, infant bonnet, red	Size 6 (34 - 36 cm)	10
4907005	nFlow CosyFit, infant bonnet, green	Size 5 (32 - 34 cm)	20
4907004	nFlow CosyFit, infant bonnet, light orange	Size 4 (30 - 32 cm)	20
4907003	nFlow CosyFit, infant bonnet, blue	Size 3 (28 - 30 cm)	20
4907002	nFlow CosyFit, infant bonnet, yellow	Size 2 (26 - 28 cm)	20
4907001	nFlow CosyFit, infant bonnet, dark grey	Size 1 (24 - 26 cm)	20
4907110	nFlow CosyFit, infant bonnet, pink	Size 0 (22 - 24 cm)	20
4907100	nFlow CosyFit, infant bonnet, grey	Size 00 (20 - 22 cm)	20
4907000	nFlow CosyFit, infant bonnet, white	Size 000 (18 - 20 cm)	20
<b>C nFlow nasal masks</b>			
		<b>Size</b>	
4716000	nFlow nasal mask	Extra large	30
4706000	nFlow nasal mask	Large	30
4705000	nFlow nasal mask	Medium	30
4704000	nFlow nasal mask	Small	30
<b>D nFlow nasal prongs</b>			
		<b>Size</b>	
4703000	nFlow nasal prong	Large	30
4702000	nFlow nasal prong	Medium	30
4701000	nFlow nasal prong	Small	30
4711000	nFlow nasal prong	Extra small	30

Fabian is a trade name of Acutronic Medical Systems.

IS1.59 BENL ▪ issue 1 04.22



Intersurgical Benelux B.V. Vluchttoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016, ISO 14001:2015 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.