


# STARMed CaStar 'R' NEXT

Hood without underarm straps for non invasive mechanical ventilation (NIV)



 **INTERSURGICAL**<sup>®</sup>  
COMPLETE RESPIRATORY SYSTEMS  
Quality, innovation and choice

[www.intersurgical.com](http://www.intersurgical.com)



@intersurgical



/intersurgical

# STARMed CaStar 'R' Next hood for non invasive ventilation (NIV)

## StarMed CaStar R Next

The CaStar R Next is a comfortable, versatile and lightweight patient interface for non invasive mechanical ventilation (NIV) designed to offer an alternative to traditional methods of NIV delivery. The CaStar R Next is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in both semi-recumbent and supine positions providing an ideal option for NIV delivery.

The CaStar R Next has been designed to improve patient-ventilator interaction, by reducing the inspiratory trigger delay and increasing the rate of airway pressurisation, while

avoiding the need for underarm straps.

The CaStar R Next may be used for the treatment of a number of respiratory conditions including:


Pre and post extubation respiratory support, COPD exacerbation, thoracic trauma, asthma, post operative hypoxaemia, cardiogenic and non cardiogenic pulmonary oedema.

**Reference:** Intensive Care Med. 2013 Spr39(4): 734-8  
*Bench comparative evaluation of a new generation and standard helmet for delivering non-invasive ventilation.*  
 Oliveri C, Costa R, Spinazzola G, Ferrone G, Longhini F, Cammarota G, Conti G, Navalesi P.

## Features and Benefits



CE certified product



**New fastening system**  
 Replaces underarm straps



**New sealed access**  
 For probes or catheters

CaStar R Next hood with anti-suffocation valve, patient access port, ear plugs and measuring tape

Code	Size	Neck Circumference	Box qty.
CP238XXL/2RW	Extra Extra Large	48-55cm	5
CP238XL/2RW	Extra Large	42-48cm	5
CP238L/2RW	Large	35-42cm	5
CP238M/2RW	Medium	28-35cm	5
CP238S/2RW	Small	21-28cm	5
CP238XS/2RW	Extra Small	≤21cm	5

Single packed. Minimum purchasing package: 5 units. Custom-packs are supplied on specific request of the customer.

IS1.24\_UK • Issue 3\_05.13



Crane House, Molly Millars Lane,  
 Wokingham, Berkshire, RG41 2RZ  
 T. +44 (0)118 9656 300  
 F. +44 (0)118 9656 356

UK  
[info@intersurgical.co.uk](mailto:info@intersurgical.co.uk)  
 Italy  
[info@intersurgical.it](mailto:info@intersurgical.it)

España  
[info@intersurgical-es.com](mailto:info@intersurgical-es.com)  
 Portugal  
[info@intersurgical.pt](mailto:info@intersurgical.pt)

Lietuva  
[info@intersurgical.lt](mailto:info@intersurgical.lt)  
 Россия  
[info@intersurgical.ru](mailto:info@intersurgical.ru)

France  
[info@intersurgical.fr](mailto:info@intersurgical.fr)  
 South Africa  
[info@intersurgical.co.za](mailto:info@intersurgical.co.za)

Japan  
[info@intersurgical.co.jp](mailto:info@intersurgical.co.jp)  
 Nederland  
[info@intersurgical.nl](mailto:info@intersurgical.nl)



Manufactured by: Intersurgical S.p.A.