

Designed to Care



nFlow™

Infant nasal CPAP system



Critical Care ▪ Breathing Systems





nFlow™

Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

A choice for infant CPAP

The Intersurgical nFlow system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.

nFlow is compatible with flow drivers designed for use with gas flip technology CPAP, providing slightly higher system pressure associated with this treatment.

The nFlow system

The nFlow system is supplied assembled with generator, heated wire limb and 10mm adjustable Superset style exhaust tube for efficient gas removal. The complete nFlow breathing system is phthalate free including the small bore tube used for monitoring and gas delivery.


A choice of patient interfaces and bonnets

Patient interfaces and bonnets are supplied separately to avoid waste. Use the size guide provided to choose the correct option for your patient.



Make an enquiry


nFlow setup video is available at www.intersurgical.nl





Choice made easy

The nFlow system includes two sizing guides making it really simple to choose the correct fitting interface and bonnet for the comfort of your patient.




12 bonnets

Available separately in twelve colour-coded sizes. Each bonnet includes three Velcro® fasteners and secure attachment points.



Patient interface

Masks and prongs are available separately in four size options. Made from soft, comfortable silicone to ensure optimal patient fit with minimal wastage.



B Bonnets



C Nasal masks



D Nasal prongs





We've made ordering easy. To make a complete system you can choose one product from group A and B, and choose a patient interface from C and/or D. At any time you can order separate items according to your needs.



Code	Description	Box Qty.
------	-------------	----------

A nFlow breathing systems

	4700000	nFlow infant CPAP heated wire breathing system, 1.6m	12
	4700001	nFlow infant CPAP heated wire breathing system, 1.8m	10
	4700016	nFlow infant CPAP heated wire breathing system with detachable monitoring line, 1.6m	10
	4708000	nFlow silencer kit	50
	4700011	nFlow replacement bonnet ties	60

B nFlow bonnets

		Size	
	4707009	nFlow, infant bonnet, navy	size 9 (40 - 42 cm) 10
	4707008	nFlow, infant bonnet, turquoise	size 8 (38 - 40 cm) 10
	4707007	nFlow, infant bonnet, orange	size 7 (36 - 38 cm) 10
	4707006	nFlow, infant bonnet, red	size 6 (34 - 36 cm) 10
	4707005	nFlow, infant bonnet, green	size 5 (32 - 34 cm) 20
	4707004	nFlow, infant bonnet, light orange	size 4 (30 - 32 cm) 20
	4707003	nFlow, infant bonnet, blue	size 3 (28 - 30 cm) 20
	4707002	nFlow, infant bonnet, yellow	size 2 (26 - 28 cm) 20
	4707001	nFlow, infant bonnet, dark grey	size 1 (24 - 26 cm) 20
	4707110	nFlow, infant bonnet, pink	size 0 (22 - 24 cm) 20
	4707100	nFlow, infant bonnet, grey	size 00 (20 - 22 cm) 20
	4707000	nFlow, infant bonnet, white	size 000 (18 - 20 cm) 20

C nFlow nasal masks

		Size	
	4716000	nFlow nasal mask	extra large 30
	4706000	nFlow nasal mask	large 30
	4705000	nFlow nasal mask	medium 30
	4704000	nFlow nasal mask	small 30

D nFlow nasal prongs

		Size	
	4703000	nFlow nasal prong	large 30
	4702000	nFlow nasal prong	medium 30
	4701000	nFlow nasal prong	small 30
	4711000	nFlow nasal prong	extra small 30

IS1.16 NL • issue 6 11.16



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal
Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia
Czech Republic ▪ Turkey ▪ South Africa ▪ China ▪ Japan
Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia



Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485

Please think before you print
Save energy and paper.
If you must print this e-mail
please print it double sided