

Explorer Endoscopy masks

For analgesic gas administration during endoscopy procedures



Explorer Endoscopy masks

Endoscopy masks for single patient use to prevent the possibility of cross infection.

The Intersurgical Explorer masks are designed to allow easy insertion of a telescope probe into the nose or mouth whilst administering an analgesic gas at the same time. There are three one-way valves positioned on the mask (with the exception of the size 1 where there is only one valve), which provide easy access from either side of the patient.

The masks are principally used for bronchoscopy and endoscopy procedures, with the analgesic gas being administered via the central 22mm port. These masks can also be used on patients undergoing a transesophageal endocardiology (T.O.E) procedure.

More information?

Call **0118 9656 376**
www.intersurgical.com

Features and benefits

Colour-coded hook ring

For ease of identification, and integration with other Intersurgical airway products.

Anatomical shape

For a precise fit.

4 sizes available


From Adult to Infant

Soft cushion




For patient comfort and air tight fit.

Self-sealing valves

To facilitate easy access.



Ordering information



Code: 1591			Box qty. 40
Infant Explorer Endoscopy Mask	size 1		
Code: 1592			Box qty. 35
Paediatric Explorer Endoscopy Mask	size 2		
Code: 1593			Box qty. 50
Small Adult Explorer Endoscopy Mask	size 3		
Code: 1594			Box qty. 40
Adult Explorer Endoscopy Mask	size 4		

IS5.18 NL • Issue 3 11.12