



Right-angled water traps

Self-sealing water traps



Critical Care ▪ Accessories



Right-angled water traps

Self-sealing water traps

Introducing Intersurgical's new right-angled water traps for the collection of condensation.

Complimenting our existing range of products, the angled water trap is now available in two sizes. Medium volume 70ml and high volume 125ml, both are ideal for use in applications where a 90° attachments may be required.

The self-sealing facility allow the water traps to be emptied without interrupting gas flow.

Both water traps feature a 22F machine side connection to a 22M breathing system side connection, and can be attached directly to a ventilator or to an anaesthetic machine.



Self-sealing facility

Allow the water traps to be emptied without interrupting flow



1916000

1915000

Code	Description	Box Qty.
1916000	Water trap, right-angled, self-sealing, medium 70ml	55
1915000	Water trap, right-angled, self-sealing, large 125ml	40

IS10.33 NL • Issue 3 11.14



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK • Ireland • France • Germany • Spain • Portugal • Italy
Netherlands • Sweden • Lithuania • Russia • Czech Republic
South Africa • China • Japan • Taiwan • Philippines
USA • Colombia • Australia



Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485

Please think before you print
Save energy and paper.
If you must print this document
please print it double sided