



## QuadraLite™

A truly unique anaesthetic face mask



**Anaesthesia** • Anaesthetic Face Masks



## QuadraLite™ anaesthetic masks

22F

### Transparent shell

For excellent visibility

### Anatomical shape

Provides a precise fit combined with lower dead space

### Soft seal

A comfortable airtight fit with flexibility for different face sizes


### Colour-coding

For easy identification of size and selection




### Unique design developed by Intersurgical

QuadraLite does not use an inflatable cushion. Innovative moulding techniques have enabled the combination of two ideally suited materials allowing the creation of a face mask with an anatomical shape and soft seal. This provides a comfortable, superior airtight fit with a crystal clear shell to enhance patient visibility.



### Four sizes for all patients

The mask is lightweight and requires only minimal pressure to achieve a seal over a wide range of face sizes and shapes. Only four sizes instead of the traditional seven provide a solution for all patients.



### The greener, safer choice

QuadraLite is part of the Intersurgical Eco Range of products containing no PVC which are designed to have a low environmental impact. Masks are phthalate free to protect the patient and a serial number for patient records and traceability is etched directly onto every mask.






Video available at [www.intersurgical.nl](http://www.intersurgical.nl)



[www.intersurgical.com/info/quadralite](http://www.intersurgical.com/info/quadralite)



## Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite™ masks
Neonatal	Preterm	<2	0	 Code: 7191000
	Term	2-4		
Infant	3 months	6	1	
	6 months	8		
	9 months	9		
Paediatric	1 year	10	2	
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small adult	6 years	20	3	 Code: 7193000
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
	12 years	39		
Medium adult	13 years	44	4	
	14 years	50		
	15 years	54		
	16 years	58		
Large adult	16+ years	60-120	5	 Code: 7194000
Extra large adult	16 + years	>120	6	



Code	Description	Size	Box Qty.
7194000	QuadraLite™, large adult, anaesthetic mask, with green cushion and hook ring, 22F (replaces traditional mask size 5–6)	4	35
7193000	QuadraLite™, medium adult, anaesthetic mask, with yellow cushion and hook ring, 22F (replaces traditional mask size 3–5)	3	35
7192000	QuadraLite™, small adult, anaesthetic mask, with white cushion and hook ring, 22F (replaces traditional mask size 2–3)	2	25
7191000	QuadraLite™, infant, anaesthetic mask, with grey cushion, 15M (replaces traditional mask size 0–1)	1	30
7089000	Anaesthetic mask hook ring for both QuadraLite and EcoMask™ with 22mm port		150

Make an enquiry

Lower environmental impact product

IS5.24 • Issue 5 06.16



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal  
 Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia  
 Czech Republic ▪ Turkey ▪ South Africa ▪ China ▪ Japan  
 Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia



**Please think before you print**  
 Save energy and paper.  
 If you must print this e-mail  
 please print it double sided

Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485