



STARMed CaStar R Up

Openable hood for non-invasive mechanical ventilation (NIV)



Critical Care ▪ Patient Interfaces ▪ **STARMed**



Quality, innovation and choice

Interact with us



www.intersurgical.com



STARMed CaStar R Up openable hood for NIV

StarMed CaStar R Up

Openable CaStar R Up is a comfortable, versatile and lightweight patient interface for non-invasive mechanical ventilation (NIV) designed to offer an alternative to traditional methods of NIV delivery. The CaStar R Up is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in semi-recumbent, supine and prone positions providing an ideal option for NIV therapy.

Advantages

- Allows rapid access to the patient during nursing procedures or in case of emergency
- One operator can easily remove the upper part of the hood and have full access to the patient
- Probes and catheters can remain in place, on the patient, connected to the lower part of the hood

Technical Features

- Manufactured from phthalate-free materials
- 22M connectors compliant to ISO standard
- Built in inflatable neck cushion for added patient comfort
- Comfortable and atraumatic underarm fastening straps
- Other fastening accessories are available if required
- Two versatile sealed access ports suitable for single or multiple lumen probes or catheters
- New design patient access port. Fitted with bi-directional anti-asphyxiation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated

[Read more evidence on www.intersurgical.com](http://www.intersurgical.com)



CaStar R Up openable hood for NIV with inflatable neck cushion, bi-directional anti-asphyxiation valve and patient access port*

| Code | Size | Neck circumference | Box Qty. |
|-------------|-------------------|--------------------|----------|
| CA022XXL/2R | Extra Extra Large | ≥50cm | 5 |
| CA022XL/2R | Extra Large | 45-52cm | 5 |
| CA022L/2R | Large | 40-47cm | 5 |
| CA022M/2R | Medium | 34-41cm | 5 |
| CA022S/2R | Small | 27-34cm | 5 |
| CA022XS/2R | Extra Small | 17-27cm | 5 |

[Find out more](#)

[Make an enquiry](#)

*Tubing not included and has to be ordered separately. Please refer to the website for the range of accessories.

IS1.39_INT ▪ Issue 6 07.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



Manufacturer:
Intersurgical S.p.A.
Mirandola (MO) Italy

Please think before you print
Save energy and paper. If you must print this information sheet please print it double sided.