



STARMed CaStar R Up

Openable hood for non-invasive mechanical ventilation (NIV)



Critical Care ▪ Patient Interfaces ▪ **STARMed**



STARMed CaStar R Up openable hood for NIV

StarMed CaStar R Up

Openable CaStar R Up is a comfortable, versatile and lightweight patient interface for non-invasive mechanical ventilation (NIV) designed to offer an alternative to traditional methods of NIV delivery. The CaStar R Up is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in semi-recumbent, supine and prone positions providing an ideal option for NIV therapy.

Advantages

- Allows rapid access to the patient during nursing procedures or in case of emergency
- One operator can easily remove the upper part of the hood and have full access to the patient
- Probes and catheters can remain in place, on the patient, connected to the lower part of the hood

Technical Features

- Manufactured from phthalate-free materials
- 22M connectors compliant to ISO standard
- Built in inflatable neck cushion for added patient comfort
- Comfortable and atraumatic underarm fastening straps
- Other fastening accessories are available if required
- Two versatile sealed access ports suitable for single or multiple lumen probes or catheters
- New design patient access port. Fitted with bi-directional anti-asphyxiation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated

[Read more evidence on www.intersurgical.nl](http://www.intersurgical.nl)



CaStar R Up openable hood for NIV with inflatable neck cushion, bi-directional anti-asphyxiation valve and patient access port*

Code	Size	Neck circumference	Box Qty.
CA022XXL/2R	Extra Extra Large	≥50cm	5
CA022XL/2R	Extra Large	45-52cm	5
CA022L/2R	Large	40-47cm	5
CA022M/2R	Medium	34-41cm	5
CA022S/2R	Small	27-34cm	5
CA022XS/2R	Extra Small	17-27cm	5



www.intersurgical.nl/info/starmed

Make an enquiry

*Tubing not included and has to be ordered separately. Please refer to the website for the range of accessories.



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



Manufacturer:
Intersurgical S.p.A.
Mirandola (MO) Italy

Please think before you print

Save energy and paper. If you must print this information sheet please print it double sided.

IS1.39_NL • Issue 5 06.20