



Humidifier electrical adaptor leads







A range of leads for electrical connection of the heated wire breathing system to the humidifier.



Critical Care ▪ Humidification Accessories

Humidifier electrical adaptor leads

In a busy ITU there are many different types of equipment and consumables to deal with. In order to streamline your active humidification breathing system requirements, Intersurgical offer a variety of electrical adaptor leads enabling you to standardise on one heated wire breathing system, dual or single, which can be used on the humidifier bases shown below.

| Code | Description | Box Qty. | Humidifier |
|---------|---|----------|--|
| 5600000 |  <p>Electrical adaptor lead for single heated wire breathing systems (clover leaf connector)</p> | 1 |  <p>PMH7000, MR850™</p> |
| 5601000 |  <p>Electrical adaptor lead for dual heated wire breathing systems (clover leaf/oval connectors)</p> | 1 | |
| 5603000 |  <p>Electrical adaptor lead for single heated wire breathing systems (clover leaf connector)</p> | 1 | <p>MR730™</p>  |
| 5604000 |  <p>Electrical adaptor lead for dual heated wire breathing systems (clover leaf/oval connectors)</p> | 1 | |

MR850 and MR730 are trademarks of Fisher & Paykel Healthcare.